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© 2002, 2011, 2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, FluiDots, and MEDRAD FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the above listed sterile disposable products can be found on page 2 of these Instructions for Use.

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Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for use: The Cascade Transfer Set with dual check valve is a sterile, single use disposable device for transferring contrast media from three contrast bottles into a Front Load Syringe. This transfer set is intended for use with MEDRAD® CT power injector applications. The MEDRAD CT power injector and Front Load Syringe must be operated in accordance with their instructions for use.

Contraindications: These devices are not intended for multiple patient-use nor for any other use for which they are not indicated.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

⚠ WARNING

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Do not use with angiographic injectors or high pressure applications greater than 300 PSI (20.7 bar). Pressures generated by Angiographic injectors may result in damage to the dual check valve, causing contrast spills which could cause personal injury and/or equipment damage.

Do not use with flow rates higher than 4 milliliters per second, or pressures above 170 PSI (11.7 bar). Pressures generated by faster flow rates may result in damage to the dual check valve. A damaged dual check valve may cease to function as intended, or could leak. Leaks and spills of contrast media could cause hazards resulting in personal injury and/or equipment damage.

Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid paths. Air embolization can cause patient injury or death. Operator vigilance and care, coupled with an appropriate procedure is essential to the avoidance of air embolism. Refer to the Injector's Operator Manual and syringe "Instructions for Use" for additional information regarding loading procedures and use MEDRAD® Fluidot indicators.

Ensure fluid path is open before injecting. Failure to do so may cause damage to the dual check valve, resulting in spills which could cause personal injury and/or equipment damage.

Spikes must be connected to contrast bottles in the sequence indicated in Step 1. Failure to do so may result in contrast spills which could cause personal injury and/or equipment damage.

⚠ WARNING

All spikes must be connected to contrast bottles before use. Failure to do so may result in contrast spills which could cause personal injury and/or equipment damage.

Operator injury could result from using a sharp spike. Use care in handling and inserting the spike into the bottle.

Use aseptic technique in product handling and making connections to avoid contamination.

The PVC tubing in this product contains DEHP, a chemical known to the State of California to cause birth defects or other reproductive harm.

Extremely high doses of phthalates such as DEHP have been reported to cause birth defects in research animals. Although limited scientific data is available on the effect of phthalates in humans, health authorities recommend that exposure to devices containing these compounds be minimized as much as possible for children and nursing or pregnant women.

DEHP is also a substance on the candidate list of Substances of Very High Concern (SVHC) according to Article 59 of REGULATION (EC) no. 1907/2006 concerning the Registration, Evaluation, Authorization and Restriction of Chemicals (REACH).

Substance Name	CAS Number	Ec Number
Bis (2-ethyl(hexyl)phthalate) (DEHP)	117-81-7	204-211-0

⚠ CAUTION

Transfer set connections must be hand tight to avoid damage to the syringe and the transfer set. Use 1/4 to 1/2 turns when connecting the Cascade Transfer Set to Front Load Syringes for MEDRAD CT injectors.

Avoid or remove kinks in the tubing by adjusting the placement of the spikes. Kinked tubing will reduce the filling rate performance of the Cascade Transfer Set.

PROCEDURE FOR USE


Figure 1

- Inspect the transfer set and package for damage. Do not use if the package is opened or damaged.
- Use aseptic technique.
- The syringe must be attached to the injector and pointed upward. The piston and syringe plunger must be in a full forward position when filling the power injector syringe.
- Use slow filling rates to reduce formation of bubbles when filling the power injector syringe.

Step 1. Connect Transfer Set to the Contrast Bottles (see Warnings and Cautions).

- Prepare the contrast bottles for spiking.
- Remove the protective cover from spike #1, (see Figure 1), and insert it into a contrast bottle.
- Remove the protective cover from spike #2, (see Figure 1), and insert it into a contrast bottle.
- Remove the protective cover from spike #3, (see Figure 1), and insert it into a contrast bottle.

Step 2. Connect Transfer Set to Front Load Syringe for CT (see Cautions).

- Remove the orange protective cap from the female luer of the transfer set.
- Remove the protective cap from the Front Load Syringe.
- Connect the transfer set to the Front Load Syringe using hand-tight connections.

Dual Check Valve Operation

- When connected to the syringe, retracting the plunger draws the contents from the spiked contrast bottles into the syringe.
 - When connected to the syringe, advancing the plunger forces the contents of the syringe through the male luer of the dual check valve.
- NOTE:** The dual check valve may make a squealing sound during operation, particularly when drawing or pushing air. This is considered normal.

Cascade Transfer Set Operation

- When connected to the Cascade Transfer Set, the three contrast bottles behave as a single reservoir of contrast media.
- As contrast media is drawn through the dual check valve during syringe filling, the contrast bottles empty in a cascading manner.
- The contrast bottles are automatically emptied in succession from bottle 3 into bottle 2, and from bottle 2 into bottle 1 (see Figure 1). In this way, bottle 3 will empty first, bottle 2 will empty second, and bottle 1 will empty last.

NOTE: Flow rates may vary between the contrast bottles and the dual check valve. Flows between bottles in the three-bottle reservoir may continue after syringe filling is completed. These flows will end when pressure within the reservoir system equalizes. This is considered normal.

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