Sterile Disposable Syringe Kit Kit de seringue jetable stérile

Instructions for Use / Mode d'emploi

Syringe Kit / Kit de seringue

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©2010-2014, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, MEDRAD Stellant, Stellant, MEDRAD FluiDots, and FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the above listed sterile disposable products can be found on page 2 of these Instructions for Use.

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ENGLISH

Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for Use: The contents of this package are intended to be used in the delivery of contrast media or saline. They are indicated for single-use on one patient only with MEDRAD® Stellant Injectors.

Contraindications: These devices are not intended for multiple patient use, drug infusion, chemotherapy, or any other use for which the device is not indicated.

Restricted Sale: Rx Only

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

Warnings

Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD® FluidDots indicators (where applicable) to reduce the chance of air embolism.

Reuse of this product may result in biological contamination, product degradation and/or product performance issues. Properly discard disposable items after use, or if there is any possibility that contamination may have occurred.

Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Bacterial contamination can occur if syringes are used to store contrast media. Use loaded syringes immediately. Discard unused loaded syringes.

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed as validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Patient or operator injury may result from contrast media leaks or tubing ruptures. Ensure that the fluid path is open; do not exceed pressures identified on the front of the package. Use of greater pressures or occlusions in the fluid path may result in leaks or ruptures.

Patient injury could result if syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged.

Use care in handling and inserting spike into contrast bottle. The spike is sharp and may cause personal injury.

Contamination may occur if the end of the spike is touched.

Do not use the transfer set on more than one container of media.

Do not use the transfer set for longer than six hours, once it has been exposed to media.

Patient or operator injury may result from damage to the transfer set valve. Do not use with injection flow rates higher than 4 mL/sec. Pressure generated by faster flow rates may cause damage to the valve resulting in leaks or ruptures.

Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten. This will help minimize leaks, disconnection, and component damage.

Refer to the Injector Operation Manual for further instructions.

Note: To minimize bubbles when loading a syringe through a transfer set, do not exceed a reverse load rate of 7 mL/sec.

Installing a Syringe

Insert the syringe until it snaps into place. The piston automatically engages the plunger, and advances it to the front of the syringe (most models). To reduce the size and quantity of air bubbles drawn into the syringe during loading, a Fluid Dispensing Device from Bayer, is recommended.

Operator vigilance and care, coupled with a set procedure, is essential to minimizing the possibility of an air embolism. Point the injector head up during loading. Point the injector head down during an injection.

To help avoid an air injection, syringes from Bayer are equipped with FluidDots indicators. FluidDots indicators should be observed as part of an arming procedure. When the FluidDots are viewed through an empty syringe, the dots appear as small narrow ellipses. When viewed through a full syringe, the dots become larger, almost round. See below:

To minimize air embolization risks, ensure that one operator is designated the responsibility of filling the syringe(s). Do not change operators during the procedure. If an operator change must occur, ensure that the new operator verifies that the fluid path is purged of air.

Loading and Priming a Syringe Using a Spike or a Quick Fill Tube (QFT)

The syringe can be loaded manually or automatically on the MEDRAD Stellant injector. See the MEDRAD Stellant Operation Manual for more detailed instructions.

MEDRAD Stellant Manual Loading:

1. Install a new syringe.
2. Install the spike or QFT onto the end of the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
3. Insert the spike or QFT into the fluid source and fill the syringe with fluid, using the load buttons.
4. Expel air.
5. Remove spike or QFT. Connect the disposable tubing set.
6. Follow the instructions in the section "Connector Tube Installation."

MEDRAD Stellant Integral Autoload with a spike or QFT:

1. Set and lock protocol.
2. Install a new syringe.
3. Install the spike or QFT onto the end of the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
4. Insert the spike or QFT into the fluid source and press Autoload.
5. Press Fill.
6. Expel air.
7. Remove spike or QFT. Connect the disposable tubing set.
8. Follow the instructions in the section "Connector Tube Installation."
9. Press Start/Hold.

NOTE: Patency check through aspiration using the injector manual knob is not possible when using a check valve on the end of the connector tube. If aspiration is important, remove the check valve from the connector tube and directly connect the connector tube to the catheter.

Removing a Syringe(s)

1. Disconnect the disposable tubing set from the vascular entry device. The disposable tubing set does not need to be disconnected from the syringe.
2. Rotate the syringe approximately 1/4 turn counter-clockwise and gently pull the syringe out of the injector head, discarding the syringe with disposable tubing set.

Note: Once the syringe is removed from the injector, the piston will automatically retract (most models).

Note: In order to remove the syringe, the last piston motion must be in the forward direction, which is typical. If you cannot remove the syringe, rotate the injector head manual knob approximately one revolution in the forward direction, then repeat Step 2.

Connector Tube Installation

1. Remove the connector tube from the package, removing the dust covers on the luer fittings.
2. Ensure all air is purged from the syringe.
3. Attach the connector tube to the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
4. Ensure that the connector luer fitting is secured to the tip of the syringe and verify that the tubing is not kinked or obstructed.

If you are using a connector tube with T-connector, attach the straight portion of the T-connector to the contrast (Syringe A) and the extension to the saline (Syringe B). If the T-connector is connected to Syringe B, priming will not completely fill the connector tubing.

5. Connect the Prime Tube to the patient end of the tubing set, then prime the tubing by pressing Prime.

6. Ensure all air is purged.

7. Rotate the injector head downward.

8. Remove the Prime Tube, connect to the patient and press Check For Air.

Transfer Set Installation

1. Connect the transfer set by inserting the media container into the holder, remove the spike cover and insert the spike.
2. Remove the protective caps from the syringe and female luer of the transfer set.
3. Install the transfer set onto the end of the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
4. Fill the syringe by retracting the syringe piston/plunger. When complete, expel air from the syringe and fluid path.

Note: The transfer set valve may make a squealing sound during operation. This is considered normal.

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