IMAXEON

INSTRUCTIONS FOR USE

150-EF-Q

Angiographic Syringe and Quick Fill Tube

Catalog No. 150-EF-Q

Imaxeon Pty Ltd
Unit 1, 38-46 South Street
Rydalmere NSW 2116
Australia

Bayer Medical Care B.V.
Avenue Céramique 27
6221 KV Maastricht
The Netherlands

Language	Page
English	4
Deutsch	5
Español	6
Français	7
Italiano	8
Nederlands	9
Dansk	10
Türkçe	11
Português	12

For use only on





IMAXEON




AVIDIA™
ANGIOGRAPHIC CONTRAST INJECTOR

© 2020 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.

Bayer, the Bayer Cross, Imaxeon, MEDRAD FluiDots, MEDRAD and FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

ENGLISH	GERMAN	FRENCH	ITALIAN	PORTUGUESE (EUR)	SPANISH (EUR)	DANISH
English	Deutsch	Français	Italiano	Português	Español	Dansk
Description	Beschreibung	Description	Descrizione	Descrição	Descripción	Beskrivelse
Consult instructions for use (ISO 15223-1, 5.4.3)	Gebrauchsanleitung lesen	Consultez le mode d'emploi.	Consultare le istruzioni per l'uso	Consulte as instruções de utilização	Consultar las instrucciones de uso	Læs brugsanvisningen
Non-Pyrogenic Fluid path (ISO 15223-1, 5.6.2 & 5.6.3)	Pyrogenfreie Infusionszuleitung	Ligne d'injection aprotégée	Linea del fluido aprotgena	Caminho de fluido aprotéxico	Circuito de líquidos aprotéxico	Ikke-pyrogen væskebane
Do not use if package is opened or damaged (ISO 15223-1, 5.2.8)	Nicht benutzen, falls Packung geöffnet oder beschädigt ist	Ne pas utiliser si l'emballage est endommagé ou ouvert.	Non usare se la confezione è aperta o danneggiata	Não utilizar se a embalagem estiver aberta ou danificada	No usar si el paquete está abierto o dañado	Må ikke bruges, hvis pakken er åbnet eller beskadiget
Single Use Only (ISO 15223-1, 5.4.2)	Nur zum Einmalgebrauch.	Usage unique strict.	Esclusivamente monouso.	Para uma única utilização	Válido para un solo uso	Kun til engangsbrug.
Date of Manufacture/ Sterilisation (ISO 15223-1, 5.1.3)	Datum der Herstellung/ Sterilisation	Date de fabrication/ stérilisation	Data di fabbricazione/ sterilizzazione	Data de fabrico/ esterilização	Fecha de fabricación/ esterilización	Fremstillings- og steriliseringsdato
Use By (ISO 15223-1, 5.1.4)	Verwendbar bis	Utiliser avant le	Da usarsi entro il	Validade	Fecha de caducidad	Skal anvendes inden
Lot Number (ISO 15223-1, 5.1.5)	Partie-Nummer	Numéro de lot	Numero di lotto	Número do lote	Número de lote	Partinummer
Catalog Number (ISO 15223-1, 5.1.6)	Katalog-Nummer	Numéro de référence	Numero di catalogo	Número do catálogo	Número de catálogo	Katalognummer
Sterilised with Ethylene Oxide (ISO 15223-1, 5.2.3)	Mit Ethylenoxid sterilisiert	Stérilisé à l'oxyde d'éthylène	Sterilizzato con ossido di etilene	Esterilizado com óxido de etileno	Esterilizado con óxido de etileno	Steriliseret med ethylenoxid
Not made with natural rubber latex (ISO 15223-1, 5.4.5, Annex B, B.2)	Nicht aus Naturkautschuklatex gefertigt	N'est pas fabriqué à partir de latex de caoutchouc naturel	Non contiene lattice di gomma naturale	Não está fabricado com látex de borracha natural	No está fabricado con látex de caucho natural	Ikke fremstillet med naturligt latexgummi
Indicates the device conforms to the requirements of the European Medical Device Directive 93/42/EEC	Gibt an, daß das Gerät den Voraussetzungen der Europäische Be-stimmungen für medizinische Geräte 93/42/EEC entspricht.	Indique que le produit est conforme aux exigences de la directive européenne 93/42/CEE sur les dispositifs médicaux.	Indica la conformità alla Direttiva Europea sulle apparecchiature mediche 93/42/CEE	Indica que o dispositivo está em conformidade com os requisitos da Diretiva Europeia 93/42/CEE relativa a dispositivos médicos.	Indica que el equipo cumple los requisitos de la Directiva Europea 93/42/CEE relativa a los productos sanitarios.	Angiver at dette produkt er i overensstemmelse med kravene i European Medical Device Directive (europæiske direktiv for medicinske instrumenter) 93/42/EEC
Warning	Warnung	Avertissement	Avvertenza	Aviso	Advertencia	Advarsel
Advises you of circumstances that could result in injury or death to the patient or the operator	Umstände aufmerksam, die zu Verletzungen oder zum Tod des Patienten oder Bedieners führen könnten.	Souligne les circonstances pouvant entraîner des blessures éventuellement mortelles, pour le patient ou l'opérateur.	Avverte l'utente di circostanze lesioni anche fatali al paziente o all'operatore.	Avisa-o de circunstâncias que podem resultar em lesões ou morte do doente ou operador.	Advertencias informan de circunstancias que podrían causar lesiones o la muerte al paciente o al operador.	Oplysninger om omstændigheder, som kan føre til personskade eller død for patienten eller operatøren.
Caution	Vorsicht	Mise en garde	Attenzione	Precaução	Precaución	Forsigtig
Advises you of circumstances that could result in damage to the device	Macht auf Umstände aufmerksam, die zu Schäden am Gerät führen könnten.	Souligne les circonstances pouvant entraîner la détérioration du dispositif.	Avverte l'utente di circostanze che potrebbero danneggiare il dispositivo.	Avisa-o de circunstâncias que podem resultar em danos ao dispositivo.	Informa de circunstancias que podrían causar daños al dispositivo.	Oplysninger om omstændigheder, som kan føre til beskadigelse af instrumentet.
	Anzahl	Quantité	Quantità	Quantidade	Cantidad	Antal
	Hersteller	Fabricant	Produttore	Fabricante	Fabricante	Producent
	Zerbrechlich! Vorsichtig handhaben.	Fragile, à manipuler avec précaution	Fragile, maneggiare con cura.	Frágil, manusear com cuidado	Frágil, manipule con cuidado	Skrobelig, håndteres forsigtig
	Trocken lassen	Protéger de l'humidité	Tenere all' asciutto	Mantem seco	Mantener seco	Holdes tørt
	Diese Seite nach oben.	Cette face vers le haut	Alto	Este lado para cima	Parte superior	Denne side op
	Temperaturbereich Primäre Kennzeichnung: Lagerung Sekundäre Kennzeichnung: Transport	Plage de température Étiquette principale : Stockage Étiquette secondaire : Transport	Intervallo di temperatura Etichetta principale: conservazione Etichetta secondaria: trasporto	Intervalo de temperaturas Etiqueta principal: Armazenamento Etiqueta secundária: Transporte	Intervalo de temperatura Etiqueta principal: almacenamiento Etiqueta secundaria: transporte	Temperaturområde Primær mærkat: Opbevaring Sekundær mærkat: Transport
	Luftfugtighedsbereich Primäre Kennzeichnung: Lagerung Sekundäre Kennzeichnung: Transport	Plage d'humidité Étiquette principale : Stockage Étiquette secondaire : Transport	Intervallo di umidità Etichetta principale: conservazione Etichetta secondaria: trasporto	Intervalo de valores de humidade Etiqueta principal: Armazenamento Etiqueta secundária: Transporte	Intervalo de humedad Etiqueta principal: almacenamiento Etiqueta secundaria: transporte	Fugtighedsområde Primær mærkat: Opbevaring Sekundær mærkat: Transport
	Nicht resterilisieren	Ne pas resteriliser	Non risterilizzare	Não reesterilizar	No reesterilizar	Må ikke re-steriliseres.
	Achtung: Bitte beachten Sie die Warn- und Vorsichtsinweise in der Gebrauchsanweisung, die jedem Karton beiliegt.	Attention : Reportez-vous aux avertissements et aux mises en garde du mode d'emploi fourni dans chaque boîte.	Attenzione: Fare riferimento alle avvertenze e precauzioni riportate nelle istruzioni per l'uso incluse in ogni confezione.	Atenção: consulte os avisos e precauções nas instruções de utilização fornecidas em cada embalagem.	Atención: consulte las advertencias y precauciones de las instrucciones de uso que vienen en cada caja.	Bemærk: Der henvises til advarsler og forsigtighedsregler i brugsanvisningen i hver kation.

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Important safety notice: These devices are intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for use: These devices are intended to be used for the delivery of contrast media. These devices are indicated for single-use only with Imaxeon injectors.

Contraindications: These devices are not intended for multiple-use, drug infusion, chemotherapy, or any other use for which the devices are not indicated.

⚠️ WARNINGS

- Air embolization can cause serious patient injury or death. Do not connect a patient to the injector until all trapped air has been cleared from the syringe, connector tubing and catheter. Operator vigilance and care, coupled with a set procedure is essential for the avoidance of air embolism. See below for use of FluiDot indicators on the pressure jacket fitted to the injector.
- Refilling increases the risk of biological contamination. Do not use disposable devices on more than one patient. Discard disposable devices immediately after use.
- Patient infection may result from the use of non-sterile components. Do not remove plunger to fill syringe. Maintain sterility of all disposable components. Check for damaged packaging or components before opening.
- Patient or operator injury may result from contrast media leaks or tubing ruptures. Ensure that the fluid path is open; do not exceed pressures identified on the front of the package. Use of greater pressures or occlusions in the fluid path may result in leaks or ruptures to the connector tube.
- Filled syringes, which are stored, can promote bacterial growth. Imaxeon syringes are designed to be filled just prior to the procedure. Discard unused filled syringes.

⚠️ CAUTIONS

- Component damage may occur if not installed properly. Assure all connections are secure; do not overtighten. This will help to minimize leaks, disconnection and component damage.

General information

The 150-EF-Q syringe is a sterile, single use device. The syringe package contains a single assembled syringe and a curved filling tube, QFT, for drawing up contrast.

The syringe and filler tube have been gas sterilised, and whilst in the unopened package, retain their sterility until the USE BY date printed on the label.

The syringe is designed for use only with the Avidia Angiographic Contrast Injector - do not substitute any other type of syringe or a syringe manufactured by any alternate supplier. The syringe has been designed and tested to withstand the pressures and flow rates achievable by the AVIDIA Angiographic Contrast Injector.

Labelling information

The label of the syringe package contains important information. This will allow you to decide if the syringe is acceptable for use or should be discarded.

PLEASE CHECK THE FOLLOWING INFORMATION ON THE SYRINGE PACKAGING BEFORE USE:

- Check for use only with the Avidia Angiographic Contrast Injector, the catalog number must be 150-EF-Q.
- Check the packaging for ANY signs of damage, if in doubt discard packaging and contents.
- Check the USE BY date - discard immediately if date is exceeded.


Use information

For comprehensive use instructions, please refer to the AVIDIA Angiographic Contrast Injector Operation Manual provided with the injector.

To remove the contents from the packaging, tear the lid away from the plastic tray. Remove the syringe and curved fill tube, being careful not to touch the luer lock fittings or the fill tube ends. Inspect the syringe for any signs of damage or visible contamination, and load the syringe into the pressure jacket located at the front of the injector head.

Refer to the Operation Manual for comprehensive filling instructions.

Note: Filled syringes are easily checked by noting the shape of the FLUIDOTS on the pressure jacket. Viewed through fluid, FLUIDOTS will appear large and rounded; viewed through air or an empty syringe, the indicators will appear small and narrow in height. See the illustration below.

Expelling air and attaching the connector tube (if used)

1. After filling the syringe, expel all air from the syringe and connector tube (if used).
2. Remove the filling device and expel air by manually advancing the injector piston so that the contrast is visible in the tip of the syringe.
3. Gently tap the pressure jacket to dislodge air bubbles and expel remaining air.
4. Attach the connector tube (if used) and advance the injector piston so that contrast flows into the connector tube and is visible at the connector tip.
5. Verify that all air is expelled before connecting to the vascular entry device.

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