User Manual

MEDRAD® Spectris Solaris EP
MR Injection System Syringe Adapters

English	1	Nederlands	5	Français	9	Deutsch	15
Italiano	19						

MEDRAD® Spectris Solaris EP MR Injection System Syringe Adapters

This manual supersedes the information provided in the MEDRAD® Spectris Solaris EP (Spectris Solaris EP) Operation Manual. Please refer to this manual for Syringe adapter information.

This manual applies to the MEDRAD *Spectris Solaris EP MR Injection System Syringe Adapters*. Read all the information contained in this manual. Understanding this information will assist you in operating the device in a safe manner.

Important Safety Notice

This device is intended to be used by medical professionals with adequate training and experience in magnetic resonance imaging (MRI) studies.

Certifications

This device is designed to comply with IEC/EN 60601-1: 1988 Including Amendment 1:1991 and Amendment 2: 1995 (electrical safety).

Intended Use

The injector system is intended for the purposes of injecting intravenous MR contrast media and common flushing solutions into the human vascular system for diagnostic studies in MRI procedures. The *Spectris Solaris EP MR Injection System Syringe Adapters* are accessories to the *Spectris Solaris EP MR Injection System*.

Restricted Sale

Federal (USA) law restricts these devices to sale by or on the order of a physician.

Trademarks

Bayer, the Bayer Cross, MEDRAD, Spectris Solaris, and MEDRAD Spectris Solaris are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

Disclaimers

External wiring and modifications disclaimers: Bayer disclaims liability for any modifications or interfaces with other equipment that are not in conformity with the specifications and information contained within this manual.

Accessory equipment connected to the *Spectris Solaris EP MR Injection System* must be certified according to IEC/EN 60601-1 standard. Furthermore, all configurations shall comply with system standard IEC/EN 60601-1-1 (2nd edition) or IEC 60601-1 (3rd edition), as applicable. Anyone who connects additional equipment to the signal input or output part configures a medical system and is therefore responsible that the system complies with the requirements of the standard IEC/EN 60601-1-1 (2nd ed.) or IEC/EN 60601-1 (3rd ed.), as applicable. To obtain on-site consulting or consulting references, contact Bayer.

© 2010, 2014, 2017, 2021 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

To provide feedback or request support, please use the contact form provided on radiology.bayer.com/contact

MEDRAD® Spectris Solaris EP (Spectris Solaris EP) MR Injection System Syringe Adaptors

Through the use of syringe adaptors available only through Bayer, the MEDRAD *Spectris Solaris EP MR Injection System* is compatible with various pre-filled syringes. The system supports pre-filled syringes available from Bayer.

Warning: Patient injury can result if the adapter is not properly engaged. Ensure that the adapter is properly engaged in the injector head. Improper engagement may cause under-volume delivery.

Warning: Do not change operating personnel during the procedure. If an operator change must occur, the new operator must verify that air has been removed from the fluid path before injecting. Ensure that one operator is delegated the responsibility of installing the adapter, syringe, purging air, and arming the system.

NOTE: To avoid under-volume injections and having to repeat procedures, ensure that the proper syringe adapter is used. Refer to the Syringe Adapter Usage table in this section.

Syringe Adapter Usage

There are two color-coded adapters (green/yellow/purple) available in the Syringe Adapter Kit from Bayer, catalog #SSMR SPG.

Adapter	Syringe	Typical Flow Rates (in ml/sec)				Typical Rate Accuracy	Typical Volume Accuracy
		18 g IV Catheter	20 g IV Catheter	22 g IV Catheter	24 g IV Catheter		
		BD pn 381144	BD pn 381134	BD pn 381123	BD pn 381112		
Yellow	Bayer 15 ml	8.0	6.4	4.6	2.8	+/- (10% + 0.02 ml/s)	+/- (1% + 0.3 ml)
	Bayer 20 ml	8.0	6.4	4.6	2.8	+/- (10% + 0.02 ml/s)	+/- (1% + 0.3 ml)
Green	Bayer 5 ml	5.0	5.0	4.6	2.8	+/- (10% + 0.02 ml/s)	+/- (1% + 0.3 ml)
	Bayer 10 ml	6.0	5.5	4.6	2.8		
Purple	Eisai 13 ml	8.5	7.0	5.0	3.5		
	Eisai 17 ml	8.5	7.0	5.0	3.5	+/- (10% + 0.02 ml/s)	+/- (1% + 0.3 ml)


Installation of a Syringe Adapter

NOTE: Always visually inspect the syringe adapters for signs of deterioration or fatigue. The adapters should be replaced at the first sign of physical deterioration.


1. Install the battery and fully reverse piston A at the injector head.

NOTE: Ensure that the syringe adapter piston is positioned at the rear limit of the adapter barrel.

2. Align the flanges on the rear of the adapter with the notches in the injector head. (The adapter is keyed to properly fit only one way.)

3. Rotate the adapter 1/4 turn clockwise until it snaps into place.
4. Remove the handle from the pre-filled syringe (if applicable), then position the tip of the pre-filled syringe in the recess at the end of the syringe adapter.

5. Snap the syringe into the adapter.
6. Ensure that the syringe is fully inserted, then snap the Adapter Ring over the syringe.
7. Carefully advance the injector head piston until it comes into contact with the adapter piston.

NOTE: Advance the piston carefully. The pre-filled syringe cannot be retracted.

8. Proceed with the attachment of syringe B, disposable connections, priming of fluid path, programming, verification that air is expelled, and arming.

NOTE: If the purple syringe adapter is attached, the touch screen will display a menu prompting the user to select the type of pre-filled syringe that is attached. The user must select the appropriate pre-filled syringe to proceed with the arming sequence. The injector will not arm until a pre-filled syringe type is selected. Reverse motion and other events that remove a "Check-For-Air" confirmation will also remove a user's pre-filled syringe selection.

After the injection is complete

1. When the “Injection Complete” window is acknowledged after an injection, the piston in syringe A will automatically retract at high reverse speed if a syringe adapter is attached and volume remaining is less than 5 ml. If volume remaining is 5 ml or greater, the piston must be retracted manually.

NOTE: Auto Retract will not activate if multi-arm is selected.

NOTE: Auto Retract can be aborted by pressing any injector head key.

2. If Auto Retract is not used, manually retract the syringe A piston.
3. Open the Adapter Ring.
4. Remove the syringe by pressing on the syringe barrel through the opening in the bottom on the syringe adapter.
5. Remove syringe B.
6. Discard disposables and syringes.

Cleaning

CAUTION: Improper or careless cleaning methods may result in equipment damage. Do not soak or immerse the adapter in water. While cleaning the adapter, avoid allowing any water to leak inside system components.

- Do not use strong industrial cleaning agents or solvents such as acetone. Gently clean the adapter using a paper towel moistened with warm water and a mild disinfectant such as antibacterial hand soap.
- Check all System Safety and Warning Labels for legibility. Ensure that the labels are not damaged or missing.

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
January 25, 2021

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

The patient data that appear in this document are fictitious protected health information (PHI) and/or actual PHI from which all personally identifiable information (PII) have been removed or otherwise anonymized. No personally identifiable information is shown.

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©2010, 2014, 2017, 2021 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer Medical Care Inc.
U.S.A.
Bayer Medical Care B.V.