Addendum to ISI2 Module Installation & User Manual

**ISI2 Module for MEDRAD® Stellant CT Injection System with
Certegra® Workstation**

&

**MEDRAD® Stellant FLEX CT Injection System with Certegra®
Workstation**

medRAD® Stellant
CT Injection System

medRAD® Stellant FLEX
CT Injection System

Languages

ENGLISH	1
ČESKY	7
DANSK	13
DEUTSCH	19
EESTI	25
ESPAÑOL	31
FRANÇAIS	37
HRVATSKI	43
INDONESIA	49
ÍSLENSKA	55
ITALIANO	61
LATVIEŠU	67
LIETUVIŲ	73
MAGYAR	79
NEDERLANDS	85
NORSK	91
POLSKI	97
PORTUGUÊS (BRASILEIRO)	105
PORTUGUÊS (EUROPEU)	111
ROMÂNĂ	117
SLOVENSKINA	125
SRPSKI	131
SLOVENSKY	137
SUOMI	143
SVENSKA	149
TÜRKÇE	155
TIẾNG VIỆT	161
ΕΛΛΗΝΙΚΑ	167
БЪЛГАРСКИ	173
ҚАЗАҚША	179
МАКЕДОНСКИ	187
РУССКИЙ	193
УКРАЇНСЬКА	201
ภาษาไทย	207
한국어	213
日本語	219
简体中文	225
繁體中文	231

ISI2 Module for MEDRAD® Stellant with Certegra® Workstation CT Injection System & MEDRAD® Stellant FLEX with Certegra® Workstation CT Injection System

Addendum to Installation & User Manual

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

1 Using ISI2 Module with MEDRAD® Stellant with Certegra® Workstation CT Injection System or MEDRAD® Stellant FLEX with Certegra® Workstation CT Injection System

This document is an addendum to the ISI2 Module Installation and User Manual and provides specific information for configuring and using the ISI2 Module with a MEDRAD® Stellant with Certegra® Workstation CT Injection System or a MEDRAD® Stellant FLEX with Certegra® Workstation CT Injection System. Refer to the ISI2 Module Installation and User Manual for additional information.

NOTE: The injection system does not control the imaging system via the ISI2 Module. The ISI notifies the scanner of the injector status which enables the scanner to synchronize the scan timing based upon when the injection started. The scanner controls the start of the scan sequence after it has received the injection start status and will not start the scan unless the scanner is in the correct state. The scanner system maintains full control of the radiation initiation.

NOTE: The scanner CANNOT override any injector operation that is considered safety critical; for example, check for air, hold during an injection or stopping an injection.

Before using the instructions contained in this document:

- ◆ All physical connections must be made between the Injection system, the ISI2 Module and the scanner. Refer to the ISI2 Module Installation and User manual for installation instructions.
- ◆ Both the scanner and injector must be configured to work together.

1.1 Configure Connectivity Settings

Once all physical connections are made between the Control Room Unit (CRU), ISI2 box and scanner, the injection system software must be configured to allow communication.

1. Go to SETUP > SYSTEM SETUP and select ISI. Refer to Table 1-1: ISI Connectivity Settings

Table 1 - 1: ISI Connectivity Settings

Value	Scanner Manufacturer
OFF	(Default) ISI2 Module not active.
Legacy ISI (700/800/900)	Existing legacy ISI module (ISI 700, ISI 800, ISI 900) in use.
ISI2:CAN,ID:1	ISI 900 connection for GE with CAN stack 2.0 (Similar to ISI900G)
ISI2:CAN,ID:2	ISI 900 connection for Siemens, Toshiba (Canon), Hitachi, United Imaging with CAN stack 2.0 (Similar to ISI900S/ ISI900T, ISI900H, ISI900U)
ISI2:CAN,ID:3	ISI 900 connection for all vendors (GE, Siemens, Toshiba (Canon), Hitachi, United Imaging with CAN stack 2.4).
ISI2:DIO	ISI 700 equivalent connection for Philips scanner (SAS).

NOTE: To operate the Workstation without using the ISI2 Module, follow step 1 and select "Off"

NOTE: There are two USB ports on the ISI2 box. If the cable is connected via the wrong port for the selected configuration, a popup displays; unplug and switch ports or change the configuration.

NOTE: The system will display a popup message if the wrong ISI Connectivity value is selected for the connected scanner or ISI module.

1.2 Connection Status Indicators

The bottom right status area on the injector screens shows an icon that reflects the current connection status.

Table 1 - 2: Connection Status Indicators

Icon	Description
None	ISI2 Configuration is disabled.
	No communication between the injector and the ISI2 module.
	The user is not able to arm the injector. The injector is configured for ISI, and it is communicating to the ISI2 module. The injector and ISI2 module are functional, but the scanner signals are in a state that will not allow arming the injector.
	The injector is configured for ISI and it is communicating to the ISI2 module. The scanner signals are in a state that will allow arming the injector.
	This icon indicates the injector and scanner are functioning independently. <ul style="list-style-type: none"> ◆ When injector symbol is gray, injector is not ready and test inject will not proceed. ◆ When icon is highlighted in yellow and flashing, Injector is armed and ready for the injection. Injector and scanner communication working.
	Test inject is programmed, test inject can be started only from the injector. <ul style="list-style-type: none"> ◆ When either or both system symbol (Injector and Scanner) is gray, that system is not ready and injection will not proceed. ◆ When injector is highlighted in yellow, injector is ready and the test inject can proceed. After the test injection is completed, one of the three arrows, as shown below, will appear.
	The procedure (injection plus scan) can only be initiated by pressing the start button on the injector. <ul style="list-style-type: none"> ◆ When either or both system symbol (Injector and Scanner) is gray, that system is not ready and injection will not proceed. ◆ When both of symbols are highlighted in yellow and flashing, both systems are ready and the injection can proceed.
	The procedure (injection plus scan) can only be initiated by pressing the start button on the scanner. <ul style="list-style-type: none"> ◆ When either or both system symbol (Injector and Scanner) is gray, that system is not ready and injection will not proceed. ◆ When both of symbols are highlighted in yellow and flashing, both systems are ready and the injection can proceed.
	The procedure (injection plus scan) can be initiated by pressing the start button on EITHER the scanner or the injector. <ul style="list-style-type: none"> ◆ When either or both system symbol (Injector and Scanner) is gray, that system is not ready and injection will not proceed. ◆ When both of symbols are highlighted in yellow and flashing, both systems are ready and the injection can proceed.

NOTE: Once the injector is armed, the blinking pattern of the lights on the injector head is different depending upon whether or not ISI2 has been enabled. The lights blink faster if ISI2 is enabled than they do if ISI2 is not enabled.

NOTE: Test inject is always initiated at the injector and can be initiated regardless of the connection status with the scanner.

1.3 Troubleshooting Tips

Table 1 - 3: Troubleshooting Tips

Problem	Condition	Potential Solution
Injector and scanner will not work together	ISI2 communication icon is not present.	Enter setup and enable ISI2 Module.
	Software displays the following icon: 	There is no communication with the ISI2 Module <ul style="list-style-type: none"> ◆ Verify the cable between the injection system and the ISI2 Module is securely fastened ◆ Reboot the injector If the problem persists, contact Bayer.
	Software displays the following icon: 	There is communication with the ISI2 Module but not with the scanner. <ul style="list-style-type: none"> ◆ Verify the scanner is configured to control the injector ◆ Verify the cable between the ISI2 Module and the Scanner is securely fastened ◆ Consider rebooting the injector and/or scanner If the problem persists, contact Bayer.
	Software displays the following icon: 	<ul style="list-style-type: none"> ◆ Verify that the scanner is configured to control the injector ◆ Consult your scanner manual
The scanner cannot program the injection protocol.	(see above)	(see above actions)
Start button does not start injector when it is armed.	Software displays the following icon: 	If the arrow between the scanner gantry and injector icon points from the scanner to the injector - the injection must be started at the scanner.
Protocol resets after an injection procedure.	N/A	This is controlled by the “Reset Protocol” option in setup. Set to “No” to retain the protocol.
Protocol does not reset after an injection procedure.	N/A	This is controlled by the “Reset Protocol” option in setup. Set to “Yes” to retain the protocol.
	Was injection stopped abnormally?	“Reset Protocol” behavior does not affect abnormal terminations.

Table 1 - 3: Troubleshooting Tips

Problem	Condition	Potential Solution
Injector loses communication with the ISI2 module.	<p>After the injector is armed, software displays the following message:</p> <p>Injector lost communication with the ISI. To continue without scanner connectivity, disable ISI in Setup/System Setup.</p>	<ul style="list-style-type: none"> ◆ Disable ISI in Setup/System Setup and proceed with the injection without scanner connectivity. ◆ Verify the cable between the injection system and the ISI2 Module is securely fastened. ◆ Reboot the injector <p>If the problem persists, contact Bayer.</p>
Injector loses communication with the scanner.	<p>After the injector is armed, software displays the following message:</p> <p>“Injector lost communication with the scanner. To continue without scanner connectivity press OK.”</p>	<ul style="list-style-type: none"> ◆ Press OK to proceed with the injection without scanner connectivity. ◆ Verify the cables between the injection system and the ISI2 Module and between ISI2 Module and the scanner is securely fastened. ◆ Reboot the injector ◆ Disable ISI2 in injector software. <p>If the problem persists, contact Bayer.</p>

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
(10) B

2021-03-08

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

The patient data that appear in this document are fictitious protected health information (PHI) and/or actual PHI from which all personally identifiable information (PII) have been removed or otherwise anonymized. No personally identifiable information is shown.

Bayer, the Bayer Cross, MEDRAD, Stellant, Stellant FLEX, MEDRAD Stellant and MEDRAD Stellant FLEX, are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2021 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

To provide feedback or request support, please use the contact form provided on radiology.bayer.com/contact

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Bayer Medical Care Inc.
U.S.A.
Bayer Medical Care B.V.