MEDRAD® MRXperion

EN Instructions for Use **BG** Инструкции за употреба **CS** Návod k použití **CT** 使用説明 **DA** Brugsanvisning **DE** Gebrauchsanleitung **EL** Οδηγίες χρήσης **ES** Instrucciones de uso **ET** Kasutusjuhend **FI** Käyttöohjeet **FR** Mode d'emploi **HR** Upute za uporabu **HU** Használati útmutató **IN** Petunjuk Penggunaan **IS** Notkunarlæðbeiningar **IT** Istruzioni per l'uso **JA** 取扱説明書 **KO** 사용 지침 **KK** Пайдалану нұсқаулығы **LT** Naudojimo instrukcija **LV** Lietošanas instrukcija **NL** Gebruiksaanwijzing **NO** Bruksanvisning **PB** Instruções de uso **PL** Instrukcja użycia **PT** Instruções de utilização **RO** Instrucţiuni de utilizare **RU** Инструкция по применению **SK** Návod na použitie **SL** Navodila za uporabo **SR** Uputstvo za upotrebu **SV** Bruksanvisning **TH** คำแนะนำการใช้งาน **TR** Kullanma Talimatı **UK** Інструкції з експлуатації **VT** Hướng dẫn Sử dụng **ZH** 使用说明

REF		   
XP 65/115VS	STERILE R  ≤ 350 psi (2410 kPa)	EN Sterile Disposable MRI Kit BG Стерилен комплект за ЯМР за еднократна употреба CS Sterilní jednorázová souprava pro MRI CT 無菌拋棄式 MRT 用具 DA Sterilt engangs MRI-sæt DE Steriles Einweg-MRT-Kit EL Στείρο αναλώσιμο κιτ MRI ES Kit de RM desechable estéril ET Steriilne ühekordne MRT-komplekt FI Steriili kertakäyttöinen MRI-sarja FR Kit IRM jetable stérile HR Sterilni jednokratni pribor za MRI HU Steril egyszer használatos c. készlet IN Kit MRI Steril Sekali Pakai IS Sæft einnota MRI-sett IT Kit RM monouso sterile JA 滅菌ディスポーザブルMRIキット KO 멸균 일회용품 MRI 키트 KK Стерильденген бір реттік MPT жинағы LT Sterilus vienkartinis MRI rinkinys LV Sterils vienreizlietojams MRI komplekts NL Steriele disposable MRI-set NO Sterilt MR-sett til engangsbruk PB Kit de IRM descartável estéril PL Jałowy, jednorazowy zestaw MR PT Kit de RMN descartável estéril RO Set IRM steril, de unică folosință RU Стерильный одноразовый набор для MPT SK Sterilná jednorázová súprava pre MR SL Sterilni komplet za enkratno uporabo MRI SR Sterilni komplet za MRI za jednokratnu primenu SV Steril MRT-sats för engångsbruk TH ชุด MRI ปรจกจากเชอขชนิดใซ้แล้วทิ้ง TR Steril Tek Kullanımlık MRI Kiti UK Стерильний одноразовий комплект для MPT VT Bộ MRI dùng một lần vô trùng ZH 一次性无菌 MRI 套件

USA ONLY: A glossary of the symbols used on the MEDRAD® MRXperion sterile disposable products can be found on page 2 of these Instructions for Use.

English	12
Deutsch	14
Français	16
Italiano	18
Español (España)	20
Svenska	22
Nederlands	24
Dansk	26
Suomi	28
Polski	30
Norsk	32
Türkçe	34
Български	36
Česky	38
Ελληνικά	40
Eesti	42
Magyar	44
한국어	46
Lietuvių	48
Latviešu	50
Português (Brasil)	52
Português (Europeu)	54
Română	56
Русский	58
Slovensky	60
繁體中文	62
Íslenska	64
Slovenščina	66
Hrvatski	68
Bahasa Indonesia	70
日本語	72
ภาษาไทย	74
Tiếng Việt	76
Sprski	78
简体中文	80
Українська	82
Қазақша	84

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[illegible]

ENGLISH

Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for Use: The contents of this package are intended to be used in the delivery of contrast media or saline. They are indicated for single use on one patient only with MEDRAD® MRXperion MR Injectors.

Contraindications: None known.

Restricted Sale: Rx only

Operating temperature: +15°C to +24°C

Operating Humidity: 20% to 90%

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

REACH compliance information can be found at www.REACH.bayer.com

Warnings

Re-using fluid containers for more than one procedure may result in biological contamination. Discard contrast and saline containers after filling syringes for a single procedure.

Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD® FluiDots indicators (where applicable) to reduce the chance of air embolism.

Reuse of this product may result in biological contamination, product degradation and/or product performance issues. Properly discard disposable items after single use, or if there is any possibility that contamination may have occurred.

Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Warnings

Bacterial contamination can occur if syringes are used as storage containers. Use loaded syringes immediately. Do not use syringes to store fluids for later use. Discard unused loaded syringes.

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Patient or operator injury may result from fluid leaks or syringe or tubing ruptures. Ensure that the fluid path is open and do not use syringes and tubing rated below 350psi (2410 kPa). An occlusion in the fluid path and/or use of syringes or tubing rated below 350psi (2410 kPa) may result in leaks or ruptures.

Patient injury could result if syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged.

Use care in handling and inserting spike into the fluid source. The spike is sharp and may cause personal injury.

Contamination may occur if the end of the spike or luer connector are touched. Do not touch the end of the spike or luer connector.

Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not over-tighten. This will help minimize leaks, disconnection, and component damage.

Cautions

Refer to the Injector Operation Manual for further instructions.

Installing New Syringes

NOTE: Use good clinical practice during the loading process to maintain sterility of the syringe kit when opening.

NOTE: Disposable shall be used by use by date as shown on each package.

1. Remove the syringes from packaging.
2. Insert the syringes by pushing them down into the receptacle until they are securely in place (an audible click is heard).
3. Remove the dust covers from the tip of each syringe.
4. Ensure that the plungers are fully advanced.


Figure 1 - 1: Installed Syringes

Filling the Syringes

The syringes on the MRXperion injector can be filled automatically or manually.

NOTE: Use the small syringe (65mL) and small spike to load contrast media into Syringe A. Use the large syringe (115mL) and large spike to load saline into Syringe B.

Filling the Syringes: Automatic

NOTE: During filling, point the injector head up. During an injection, point the injector head down.

NOTE: Auto Fill of Syringe A and Syringe B cannot be started simultaneously.

1. Set and lock protocol with the desired volumes.
2. Install new syringes.
3. Remove the small spike from the packaging. Remove the dust cover from the luer end of the spike. Install the small spike onto the contrast syringe (Syringe A). Do not install with excessive force. Remove the cover

from the tip of the spike. Insert the spike into the fluid source.

4. Press the Fill A button twice to begin Auto Fill of Syringe A.

5. Ensure the fluid path is free from excess air. Remove the spike.
6. Remove the large spike from the packaging. Remove the dust cover from the luer end of the spike. Install the large spike onto the saline syringe (Syringe B). Do not install with excessive force. Remove the cover from the tip of the spike. Insert the spike into the fluid source.

7. Press the Fill B button twice to begin Auto Fill of Syringe B.

8. Ensure the fluid path is free from excess air. Remove the spike. Discard the fluid source containers and spike.

NOTE: In steps 3 and 6, if using the Female/Female Adapter (FFA) to load fluid from a pre-filled syringe, attach the FFA onto the tip of the syringe instead of a spike. (The FFA may be purchased separately. The catalog number for this item is FFA 50.)

9. Connect the disposable tubing following the instructions in the **Connector Tube Installation** section.

Filling the Syringes: Manual

NOTE: During loading, point the injector head up. During an injection, point the injector head down.

1. Install new syringes.
2. Remove the small spike from the packaging. Remove the dust cover from the luer end of the spike. Install the small spike onto the contrast syringe (Syringe A). Do not install with excessive force. Remove the cover from the tip of the spike. Insert the spike into the fluid source.
3. Fill Syringe A by pressing the Enable Piston Control button to enable the piston controls and use the reverse piston control to fill the syringe with the desired amount of fluid. Alternatively, use the manual knob on Side A of the injector.
4. Ensure the fluid path is free from excess air. Remove the spike.
5. Remove the large spike from the packaging. Remove the dust cover from the luer end of the spike. Install the large spike onto the saline syringe (Syringe B).


Do not install with excessive force. Remove the cover from the tip of the spike. Insert the spike into the fluid source.

6. Fill Syringe B by pressing the Enable Piston Control button to enable the piston controls and use the reverse piston control to fill the syringe with the desired amount of fluid. Alternatively, use the manual knob on Side B of the injector.
7. Ensure the fluid path is free from excess air. Remove the spike. Discard the fluid source containers and spike.

NOTE: In steps 2 and 5, if using the Female/Female Adapter (FFA) to load fluid from a pre-filled syringe, attach the FFA onto the tip of the syringe instead of a spike. (The FFA may be purchased separately. The catalog number for this item is FFA 50.)

8. Connect the disposable tubing following the instructions in the **Connector Tube Installation** section.

To help avoid an air embolism, syringes from Bayer are equipped with Fluidots indicators. Fluidots indicators should be observed as part of an arming procedure. When the Fluidots are viewed through an empty syringe, the dots appear as small narrow ellipses. When viewed through a fluid-filled syringe, the dots become larger (almost round).


Operator vigilance and care, coupled with a set procedure, is essential to minimizing the possibility of air embolism. Ensure that one operator is designated the responsibility of filling the syringe(s). Do not change operators during the procedure. If an operator change must occur, ensure that the new operator verifies that the fluid path is purged of air.

Connector Tube Installation

After loading and priming the syringes, install the connector tube.

1. Remove the connector tube from the package. Remove the dust covers on the luer fittings.
2. Ensure that all air is purged from the syringes.
3. Securely attach the connector tube to the syringes as shown. Do not install with excessive force.

4. Ensure that the connector luer fitting is secured to the tips of the syringes and verify that the tubing is not kinked or obstructed.

5. Prime the tubing by pressing the Prime button.
6. Ensure that all air is purged.
7. Rotate the injector head downward.

8. Examine the fluid path for air. Confirm that the fluid path is clear of excess air by pressing the Check for Air Confirmation Button. Connect to the patient.

Removing the Syringes

1. Disconnect the disposable tubing set from the vascular entry device. The disposable tubing set does not need to be disconnected from the syringes.
2. Rotate the syringes and gently pull the syringes out of the injector head. Discard the syringes and disposable tubing set.

NOTE: Once the syringes are removed from the injector head, the piston will automatically retract.

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Bayer Medical Care Inc.
U.S.A.
Bayer Medical Care B.V.
Bayer (Schweiz) AG