MEDRAD® MRXperion

EN Instructions for Use **CN** 使用说明 **DA** Brugsanvisning **DE** Gebrauchsanleitung **ES** Instrucciones de uso **FR** Mode d'emploi **IT** Istruzioni per l'uso **JA** 取扱説明書 **KK** Пайдалану нұсқаулығы **KO** 사용 지침
NL Gebruiksaanwijzing **NO** Bruksanvisning **RU** Инструкция по применению **SV** Bruksanvisning

REF		
XP 115 VS	STERILE R ≤ 350 psi (2410 kPa)	EN Sterile Disposable Syringe Kit CN 一次性无菌针筒套件 DA Sterilt engangssprøjtesæt DE Steriles Einwegspritzenkit ES Kit de jeringa desechable estéril FR Kit de seringue jetable stérile IT Kit di siringhe monouso sterili JA 滅菌ディスポーザブルシリンジキット KK Стерильденген, бір реттік шприц жинағы KO 멸균 일회용품 주사기 키트 NL Steriele disposable spuitset NO Sterilt engangssprøytesett RU Комплект стерильного одноразового шприца SV Steril engångsprutsats

USA ONLY: A glossary of the symbols used on this sterile disposable product can be found on page 2 of these Instructions for Use.

English	6
Nederlands	8
Français	10
Deutsch	12
Italiano	14
日本語	16
한국어	18
Español (España)	20
Dansk	22
Norsk	24
Svenska	26
简体中文	28
Қазақша	30
РУССКИЙ	32

ENGLISH

Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for Use: The contents of this package are indicated to be used in the delivery of contrast media or saline. They are indicated for single use on one patient only with MEDRAD® MRXperion MR Injectors.

Contraindications: None known.

Restricted Sale: Rx only.

Operating temperature: +15°C to +24°C

Operating Humidity: 20% to 90%

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

REACH compliance information can be found at www.REACH.bayer.com

Warnings

Re-using fluid containers for more than one procedure may result in biological contamination. Discard contrast and saline containers after using syringes for a single procedure.

Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of FluiDots indicators (where applicable) to reduce the chance of air embolism.

Reuse of this product may result in biological contamination, product degradation, and/or product performance issues. Properly discard disposable items after single use or if there is any possibility that contamination may have occurred.

Syringe sterility will be compromised and patient infection may result if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Bacterial contamination can occur if syringes are used as storage containers. Use loaded syringes immediately. Do not use syringes to store fluids for later use. Discard unused loaded syringes.

Warnings

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess, or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices poses risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection, as the device has not been validated to be cleaned or re-sterilized.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged or if damaged components are used. Visually inspect contents and package before each use.

Patient or operator injury may result from fluid leaks or syringe or tubing ruptures. Ensure that the fluid path is open and use disposables that are rated to the set pressure limit. An occlusion in the fluid path and/or use of syringes or tubing rated below the set pressure limit may result in leaks or ruptures.

Patient injury could result if the syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged.

Use care in handling and inserting the spike into the fluid source. The spike is sharp and may cause personal injury.

Contamination may occur if the end of the spike or luer connector are touched. Do not touch the end of the spike or luer connector.

Follow aseptic technique when handling syringe and connector tubing components. Specifically, maintain sterility of all syringe and connector tube components.

Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not over-tighten. This will help minimize leaks, disconnection, and component damage.

Refer to the injector operation manual for further instructions.

Installing New Saline Syringe

NOTE: Use good clinical practice during the loading process to maintain sterility of the syringe kit when opening.

NOTE: Disposable shall be used by the use-by date shown on each package.

1. Remove the syringe from packaging.
2. Insert the syringe by pushing it down into the injector head until it is securely in place (an audible click is heard).
3. Ensure that the plunger is fully advanced.


Figure 1 - 1: Installed Syringe

Filling the Syringe

The syringe on the MRXperion injector can be filled automatically or manually.

NOTE: During filling, point the injector head up. During an injection, point the injector head down.

Filling the Syringe: Automatic


1. Set and lock the protocol with the desired volumes.
2. Remove the dust cover from the tip of the syringe.
3. Remove the spike from the packaging. Remove the dust cover from the luer end of the spike.
4. Install the spike onto the saline syringe (Syringe B). Do not install with excessive force.
5. Remove the dust cover from the tip of the spike. Insert the spike into the fluid source.
6. Press the B Button twice to begin Auto Fill of Syringe B.
7. Ensure the syringe is free from air. Wait 3-5 seconds after Auto Fill is complete, then remove the spike and saline source. Discard the fluid source containers and spike.

B**Filling the Syringe: Manual**

1. Remove the dust cover from the tip of the syringe.
2. Remove the spike from the packaging. Remove the dust cover from the luer end of the spike.
3. Install the spike onto the saline syringe (Syringe B). Do not install with excessive force.

4. Remove the cover from the tip of the spike. Insert the spike into the fluid source.
5. Fill Syringe B by pressing the Enable Piston Control button to enable the piston controls. Use the reverse piston control to fill the syringe with the desired amount of fluid. Alternatively, use the manual knob on Side B of the injector. 
6. Ensure the syringe is free from air. Remove the spike. Discard the fluid source containers and spike.


To help avoid an air embolism, syringes from Bayer are equipped with MEDRAD® FluiDots indicators. FluiDots indicators should be observed as part of an arming procedure. When the FluiDots are viewed through an empty syringe, the dots appear as small narrow ellipses. When viewed through a fluid-filled syringe, the dots become larger (almost round.)

Installing the Connector Tube

After loading and priming the saline syringe, install the connector tube.

1. Remove the connector tube from the package. Remove the dust covers on the luer fittings
2. Securely attach the straight portion of the T-connector tube to the prefilled contrast syringe.
3. Insert the prefilled contrast syringe (PFS) into the prefilled syringe adapter (PFA) and close the orange locking ring.
4. Insert the PFA by pushing it down into the injector head until it is securely in place (an audible click is heard).
5. Attach the extension end of the connector tubing to the saline syringe.
6. Advance the adapter piston by pressing the Enable Piston Control button and using the forward piston control.
7. Ensure that the luer connections are secure. Verify that the tubing is not kinked or obstructed.

Priming

1. Prime the tubing to the T-connector with contrast using the Enable Piston Control button and arrows or the manual knob.

2. Prime the tubing with saline by pressing the Prime button.

3. Examine the syringes and the tubing to determine the fluid path is free from air. Press the I Checked For Air button to confirm that the syringes and tubing have been inspected for the presence of air.

NOTE: If required, turn the manual knobs or repeat the priming steps to advance fluid and remove any remaining air.

Preparing to Inject

1. Rotate the injector head downward until it stops and the Volume Indicators change direction.
2. Connect to the patient.
3. Arm the injector.
4. Initiate the injection.

Removing the Disposables

1. Disconnect the connector tubing set from the patient vascular entry device. The connector tubing set does not need to be disconnected from the syringes.
2. Press the A Button once to retract the adapter piston. 
3. Rotate the saline syringe and gently pull the syringe out of the injector head. Remove the prefilled syringe. Discard the syringes and connector tubing set.

NOTE: Once the saline syringe is removed from the injector head, the piston will automatically retract.

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Bayer Medical Care Inc.
U.S.A.
Bayer Medical Care B.V.
Bayer (Schweiz) AG